


HeartMate II LVAS



Destination Therapy Trials





Figure 1. Survival Rates in Two Trials of Left Ventricular Assist Devices (LVADs) as Destination Therapy.

The curves labeled 2009 are those reported by Slaughter and colleagues in this issue of the *Journal*²; those labeled 2001 were reported for the REMATCH trial.¹


Comparison of HM I (XVE) and HM II



	HM I	HM II
Weight (gm)	1250	280
Volume (ml)	450	63
Noise	Audible	Silent
Moving parts	Many	One
Maximal flow (l/min)*	10	10
Clinical Durability (yr)	1.5	Est. > 5

* at mean pressure=100 mm Hg

HM II with controller and batteries



HM II Pump External View

